

## Product Datasheet

### Human SUMO2 Protein (orb428678)

<b>Description</b>	Recombinant of human SUMO2 protein
<b>Preservatives</b>	The SUMO2 (1mg/ml) containing 20mM Tris-HCl buffer (pH 8.0).
<b>Form/Appearance</b>	Sterile filtered colorless solution.
<b>Storage</b>	Stability: Can be stored at +4C for 1 week. For long term storage , below -2°C.Please prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	Recombinant & Natural Proteins
<b>Protein Sequence</b>	MADEKPKEGVKTENNDHINLKVAGQDGSVVQFKIKRHTPLSKLMKAYCERQGLS MRQIRFRFDGQPINETDTPAQLEMEDEDTIDVFQQQTGG
<b>Purity</b>	Greater than 95.0% as determined by SDS-PAGE.
<b>Source</b>	Escherichia Coli
<b>Expiration Date</b>	12 months from date of receipt.