

## Product Datasheet

### Human TARBP2 Protein (orb428734)

<b>Description</b>	Recombinant of human TARBP2 protein
<b>Preservatives</b>	TARBP2 protein solution (1mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 10% glycerol and 0.4M Urea.
<b>Form/Appearance</b>	Sterile Filtered colorless solution.
<b>Storage</b>	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	Recombinant & Natural Proteins
<b>Protein Sequence</b>	MLAANPGKTP ISLLQEYGTR IGKTPVYDLL KAEGQAHQPN FTFRVTVGDT SCTGQGPSKK AAKHKAAEVA LKHLKGG SML EPALEDSSSF SPLDSSLPED IPVFTAAAAA TPVPSVVLTR SPPMELQPPV SPQQSECNPV GALQELVVQK GWRLPEYTVT QESGPAHRKE FTMTCRVERF IEIGSGTSKK LAKRNAAAKM LLRVHTVPLD ARDGNEVEPD DDHFSIGVGS RLDGLRNRGP GCTWDSL RNS VGEKILSLRS CSLGSLGALG PACCRVLSL SEL SEEQAFHVS Y LDIEELSLSG LCQCLVELST QPATVCHGSA TTREAARGEA ARRALQYLKI MAGSK
<b>Purity</b>	Greater than 85.0% as determined by SDS-PAGE.
<b>Source</b>	Escherichia Coli
<b>Expiration Date</b>	12 months from date of receipt.