

Product Datasheet

Human TGIF2LX Protein (orb428820)

Description	Recombinant of human TGIF2LX protein
Preservatives	The TGIF2LX solution (0.5mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.2M NaCl, 50% glycerol and 1mM DTT.
Form/Appearance	Sterile Filtered clear solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMEAADG PAETQSPVEK DSPAKTQSPA QDTSIMSRNN ADTGRVLALP EHKKKRKGNL PAESVKILRD WMYKHRFKAY PSEEEKQMLS EKTNLSLLQI SNWFINARRR ILPDMLQRR NDPIIGHKTG KDAHATHLQS TEASVPAKSG PSGPDNVQSL PLWPLPKGQM SREKQDPES APSQKLTGIA QPKKKVKVSV TSPSSPELVS PEEHADFFSF LLLVDAAVQR AAELELEKKQ EPNP
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.