

Product Datasheet

TNF a (His) Protein (orb428887)

Description	Recombinant of TNF a (His) protein
Preservatives	TNF a protein solution (1mg/ml) containing Phosphate Buffer Saline (pH7.4), 10% glycerol and 1mM DTT.
Form/Appearance	Sterile Filtered clear solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMLRSSS QNSSDKPVAH VVANHQAEQ LEWLSQRANA LLANGMDLKD NQLVVPADGL YLIYSQVLFK GQGCPDYVLL THTVSRFAIS YQEKVSLLSA IKSPCKDTP EGAELKPWYE PMYLGGVFQL EKGDLLSAEV NLPKYLDITE SGQVYFGVIA L
Purity	Greater than 90.0% as determined by analysis by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.