

Product Datasheet

Human UCHL3 Protein (orb429104)

Description	Recombinant of human UCHL3 protein. Ubiquitin carboxyl-terminal hydrolase isozyme L3 belongs to a ge...
Immunogen	non-glycosylated polypeptide chain containing 250 amino acids (1-230aa)
Preservatives	20mM Tris-HCl buffer (pH8.0), 1mM DTT, and 10% glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at 4°C if the entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage, it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid freeze thaw cycles.
Tag	20 amino acid His-tag at N-terminus
Note	For research use only
Application notes	Enzymes
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MEGQRWLPLE ANPEVTNQFL KQLGLHPNWQ FVDVYGMDPE LLSMVPRPVC AVLLLPITE KYEVFRTEEE EKIKSQGQDV TSSVYFMKQT ISNACGTIGL IHAIANNKDK MHFESGSTLK KFLEESVSMS PEERARYLEN YDAIRVTHET SAHEGQTEAP SIDEKVDLHF IALVHVDGHL YELDGRKPPF INHGETSDET LLEDAIEVCK KFMERDPDEL RFNAIALSAA
Purity	> 95% as determined by SDS-PAGE
Source	Escherichia Coli
MW	28.3kDa
Expiration Date	6 months from date of receipt.