

Product Datasheet

Feline UMOD Protein (orb429124)

Description	Feline UMOD protein
Preservatives	The UMOD protein was lyophilized from 0.4µm filtered solution at a concentration of 0.5mg/ml containing deionized water.
Form/Appearance	Filtered White lyophilized (freeze-dried) powder.
Storage	Stability: Store lyophilized protein at -20°C. Aliquot the product after reconstitution to avoid repeated freezing/thawing cycles. Reconstituted protein can be stored at 4°C for a limited period of time; it does not show any change after two weeks at 4°C
Note	For research use only
Application notes	Enzymes
Purity	Greater than 95% as determined by SDS-PAGE.
Source	Feline Urine
Solubility (25°C)	Add deionized water to prepare a working stock solution of approximately 0.5mg/mL and let the lyophilized pellet dissolve completely. Product is not sterile! Please filter the product by an appropriate sterile filter before using it in the cell culture.
Expiration Date	12 months from date of receipt.