

## Product Datasheet

### Human YWHAB (His) Protein (orb429271)

<b>Description</b>	Recombinant of human YWHAB (His) protein
<b>Preservatives</b>	The YWHAB solution (1mg/1ml) contains phosphate buffered Saline (pH7.4) and 10% glycerol.
<b>Form/Appearance</b>	Sterile Filtered colorless solution.
<b>Storage</b>	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	Protein Kinases
<b>Protein Sequence</b>	MGSSHHHHHH SGLVPRGSH MGSHMTMDKS ELVQKAKLAE QAERYDDMAA AMKAVTEQGH ELSNEERNLL SVAYKNVVG A RRSSWRVISS IEQKTERNEK KQQMGKEYRE KIEAELQDIC NDVLELLDKY LIPNATQPES KVFYLMKMGD YFRYLSEVAS GDNKQTTVSN SQQAYQEAFE ISKKEMQPTH PIRLGLALNF SVFYEILNS PEKACSLAKT AFDEAIAELD TLNEESYKDS TLIMQLLRDN LTLWTSENQG DEGDAGEGEN
<b>Purity</b>	Greater than 90.0% as determined by SDS-PAGE.
<b>Source</b>	Escherichia Coli
<b>Expiration Date</b>	12 months from date of receipt.