

Product Datasheet

Human STX1A Protein (orb429481)

Description	Recombinant of human Syntaxin 1A protein
Preservatives	20mM Tris-HCl pH7.5, 10% glycerol, and 1mM DTT
Form/Appearance	Sterile filtered colorless solution
Storage	Store at 4°C if the entire vial will be used within 2-4 weeks. Store frozen at -20°C for longer periods of time. For long term storage, it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Tag	N-terminal His-Tag
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MKDRTQELRT AKDSDDDDDDV AVTVDRDRFM DEFFEQVEEI RGFIDKIAEN VEEVKRKHSA ILASPNPDEK TKEELEELMS DIKKTANKVR SKLKSIEQSI EQEEGLNRSS ADLRIRKTOH STLSRKFEV MSEYNATQSD YRERCKGRIQ RQLEITGRIT TSEELEDMLE SGNPAIFASG IIMDSSISKQ ALSEIETRHS EIIKLENSIR ELHDMFMDMA MLVESQ
Source	E. coli
MW	26.1kDa
Expiration Date	6 months from date of receipt.