

Product Datasheet

TdT Assay Kit (orb429551)

Description	TdT Assay Kit
Reactivity	Human
Conjugation	Other
Tested Applications	ICC, IF
Target	TdT
Preservatives	All components contain 0.1% Sodium Azide (NaN ₃)Components 2, 2C and 3 contain 1% Bovine Serum AlbuminComponent M contains 80% glycerol. Phosphate buffered saline
Form/Appearance	Kit
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Expiration Date	12 months from date of receipt.