

Product Datasheet

Mouse IgG1/IgG2a Negative control (FITC/RPE) (orb429577)

Description	Mouse IgG1/IgG2a Negative control (FITC/RPE)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	FITC/RPE
Tested Applications	FC
Target	MOUSE IgG1/IgG2a NEGATIVE CONTROL
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Dual Colour combination consisting of FITC conjugated and RPE conjugated monoclonal antibodies mixed in optimal ratio - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Cocktail
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.