

Product Datasheet

Mouse IgG2a Negative control (FITC) (orb429591)

Description	Mouse IgG2a Negative control (FITC)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Human T lymphocytes.
Target	MOUSE IgG2a NEGATIVE CONTROL
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.