

Product Datasheet

Rat IgG2a Negative control (orb429600)

Description	Rat IgG2a Negative control
Species/Host	Rat
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human lymphocytes.
Target	RAT IgG2a NEGATIVE CONTROL
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.