

Product Datasheet

Rat IgG2a Negative control (APC) (orb429621)

Description	Rat IgG2a Negative control (APC)
Species/Host	Rat
Reactivity	Rat
Conjugation	APC
Tested Applications	FC
Immunogen	Trinitrophenol
Target	RAT IgG2a NEGATIVE CONTROL
Preservatives	0.09% Sodium Azide (NaN ₃)1% Bovine Serum Albumin5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Allophycocyanin (APC) - Lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.