

Product Datasheet

Rat IgG2b Negative control (FITC) (orb429631)

Description	Rat IgG2b Negative control (FITC)
Species/Host	Rat
Reactivity	Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	KLH
Target	RAT IgG2b NEGATIVE CONTROL
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.