

Product Datasheet

Mouse IgG1 Negative control (APC) (orb429644)

Description	Mouse IgG1 Negative control (APC)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	APC
Tested Applications	FC
Target	MOUSE IgG1 NEGATIVE CONTROL
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to APC-lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.