

Product Datasheet

Mouse IgG2a Negative control (Biotin) (orb429651)

Description	Mouse IgG2a Negative control (Biotin)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC
Immunogen	Activated rat T-helper cells.
Target	MOUSE IgG2a NEGATIVE CONTROL
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.