

Product Datasheet

Mouse IgG2a Negative control (RPE) (orb429654)

Description	Mouse IgG2a Negative control (RPE)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	RPE
Tested Applications	FC
Immunogen	Activated rat T-helper cells.
Target	MOUSE IgG2a NEGATIVE CONTROL
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.