

Product Datasheet

ALPHA MELANOCYTE STIMULATING HORMONE Peptide (orb429667)

Description	ALPHA MELANOCYTE STIMULATING HORMONE Peptide
Tested Applications	ELISA
Target	ALPHA MELANOCYTE STIMULATING HORMONE
Preservatives	None
Form/Appearance	Purified synthetic peptide - liquid
Concentration	Total protein concentration 1 mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Purity	Purified
Expiration Date	12 months from date of receipt.