

Product Datasheet

CLOSTRIDIUM DIFFICILE TOXIN A antibody (orb429725)

Description	Rabbit polyclonal antibody to CLOSTRIDIUM DIFFICILE TOXIN A
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Synthetic peptide corresponding to the amino acids 2123 - 2132 of the native molecule conjugated to Keyhole limpet Haemocyanin. Peptide sequence: CTIDGKKYYFN.
Target	CLOSTRIDIUM DIFFICILE TOXIN A
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.