

## Product Datasheet

### CD121a antibody (orb429734)

|                            |                                                                                                                   |
|----------------------------|-------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to CD121a                                                                               |
| <b>Species/Host</b>        | Mouse                                                                                                             |
| <b>Reactivity</b>          | Human                                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                                      |
| <b>Tested Applications</b> | ELISA, WB                                                                                                         |
| <b>Immunogen</b>           | Recombinant peptide corresponding to the extracellular domain of human CD121a expressed in insect cell line Sf21. |
| <b>Target</b>              | CD121a                                                                                                            |
| <b>Preservatives</b>       | 5% Trehalose. Phosphate buffered saline                                                                           |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                                                                             |
| <b>Concentration</b>       | 0.5 mg/ml                                                                                                         |
| <b>Storage</b>             | -20°C only (ship +4°C)                                                                                            |
| <b>Note</b>                | For research use only                                                                                             |
| <b>Isotype</b>             | IgG1                                                                                                              |
| <b>Clonality</b>           | Monoclonal                                                                                                        |
| <b>Purity</b>              | Purified                                                                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                   |