

Product Datasheet

ALDOSTERONE antibody (orb429739)

Description	Mouse monoclonal antibody to ALDOSTERONE
Species/Host	Mouse
Reactivity	Aves
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Aldosterone conjugated to bovine serum albumin.
Target	ALDOSTERONE
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.