

## Product Datasheet

### IL-4 RECEPTOR antibody (orb429746)

|                            |                                                                                                                                    |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to IL-4 RECEPTOR                                                                                         |
| <b>Species/Host</b>        | Mouse                                                                                                                              |
| <b>Reactivity</b>          | Human                                                                                                                              |
| <b>Conjugation</b>         | Unconjugated                                                                                                                       |
| <b>Tested Applications</b> | ELISA, FA, FC, IHC-P, WB                                                                                                           |
| <b>Immunogen</b>           | Spodoptera frugiperda lepidopteran ovarian cell line Sf21-derived, recombinant human interleukin 4 receptor, extracellular domain. |
| <b>Target</b>              | IL-4 RECEPTOR                                                                                                                      |
| <b>Preservatives</b>       | 0.5% Trehalose. Phosphate buffered saline                                                                                          |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                                                                                              |
| <b>Concentration</b>       | 1.0 mg/ml                                                                                                                          |
| <b>Storage</b>             | -20°C only (ship +4°C)                                                                                                             |
| <b>Note</b>                | For research use only                                                                                                              |
| <b>Isotype</b>             | IgG2a                                                                                                                              |
| <b>Clonality</b>           | Monoclonal                                                                                                                         |
| <b>Purity</b>              | Purified                                                                                                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                    |