

Product Datasheet

CORTICOSTERONE antibody (orb429756)

Description	Rabbit polyclonal antibody to CORTICOSTERONE
Species/Host	Rabbit
Reactivity	Aves
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Corticosterone 3 -BSA conjugate.
Target	CORTICOSTERONE
Preservatives	0.09% Sodium Azide (NaN ₃)
Form/Appearance	Phosphate buffered saline. Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Purified IgG prepared by affinity chromatography on Protein G.
Dilution Range	ELISA: 1/500-10000
Expiration Date	12 months from date of receipt.