

Product Datasheet

ACTH antibody (orb429770)

Description	Mouse monoclonal antibody to ACTH
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	Synthetic peptide corresponding to the amino acids 1 - 24 (SYSMEHFRWGKPVGKKRRPVKVYP) of the human ACTH conjugated to Keyhole Limpet Haemocyanin.
Target	ACTH
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.