

Product Datasheet

ACETYLCHOLINESTERASE antibody (orb429774)

Description	Rabbit polyclonal antibody to ACETYLCHOLINESTERASE
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Highly purified acetylcholine esterase from Electrophorus electricus. EC no. 3.1.1.7
Target	ACETYLCHOLINESTERASE
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.