

Product Datasheet

C3b antibody (orb429790)

Description	Mouse monoclonal antibody to C3b
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH-Fr, WB
Immunogen	Human C3b.
Target	C3b
Preservatives	0.09% Sodium Azide, 0.5% Bovine Serum Albumin. These concentrations relate to the reconstituted product.. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	50 ug/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.