

## Product Datasheet

### C3a antibody (orb429791)

|                            |                                                                           |
|----------------------------|---------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to C3a                                          |
| <b>Species/Host</b>        | Mouse                                                                     |
| <b>Reactivity</b>          | Human                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                              |
| <b>Tested Applications</b> | ELISA, IH-Fr, WB                                                          |
| <b>Immunogen</b>           | Human complement component C3.                                            |
| <b>Target</b>              | C3a                                                                       |
| <b>Preservatives</b>       | 0.09% Sodium Azide, 0.5% Bovine Serum Albumin.. Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified IgG - lyophilized                                                |
| <b>Concentration</b>       | 50 ug/ml                                                                  |
| <b>Storage</b>             | +4°C, -20°C if preferred                                                  |
| <b>Note</b>                | For research use only                                                     |
| <b>Isotype</b>             | IgG2b                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                |
| <b>Purity</b>              | Purified                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                           |