

Product Datasheet

ATRIAL NATRIURETIC PEPTIDE antibody (orb429792)

Description	Mouse monoclonal antibody to ATRIAL NATRIURETIC PEPTIDE
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, RIA
Immunogen	Synthetic peptide corresponding to part of the native molecule conjugated to bovine serum albumin.
Target	ATRIAL NATRIURETIC PEPTIDE
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.