

Product Datasheet

IgG antibody (orb429793)

Description	Donkey polyclonal antibody to IgG
Species/Host	Donkey
Reactivity	Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, IP
Immunogen	Sheep IgG whole molecule.
Target	IgG
Preservatives	0.01% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid (sterile filtered)
Concentration	2.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(H/L)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.