

Product Datasheet

ALBUMIN antibody (FITC) (orb429800)

Description	Sheep polyclonal antibody against ALBUMIN conjugated to FITC
Species/Host	Sheep
Reactivity	Rat
Conjugation	FITC
Tested Applications	ELISA, ICC, IF
Immunogen	Native albumin from rat serum
Target	ALBUMIN
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.