

## Product Datasheet

### IgG antibody (orb429817)

|                            |                                                                   |
|----------------------------|-------------------------------------------------------------------|
| <b>Description</b>         | Goat polyclonal antibody to IgG                                   |
| <b>Species/Host</b>        | Goat                                                              |
| <b>Reactivity</b>          | Mouse                                                             |
| <b>Conjugation</b>         | Unconjugated                                                      |
| <b>Tested Applications</b> | ELISA                                                             |
| <b>Immunogen</b>           | Mouse IgG whole molecule.                                         |
| <b>Target</b>              | IgG                                                               |
| <b>Preservatives</b>       | 0.01% Sodium Azide (NaN <sub>3</sub> ). Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified Ig - liquid                                              |
| <b>Concentration</b>       | 1.0 mg/ml                                                         |
| <b>Storage</b>             | +4°C, -20°C if preferred                                          |
| <b>Note</b>                | For research use only                                             |
| <b>Application notes</b>   | (H/L)                                                             |
| <b>Isotype</b>             | Polyclonal IgG                                                    |
| <b>Clonality</b>           | Polyclonal                                                        |
| <b>Purity</b>              | Purified                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                   |