

Product Datasheet

STREPTOCOCCUS A antibody (orb429822)

Description	Rabbit polyclonal antibody to STREPTOCOCCUS A
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Streptococci, group A
Target	STREPTOCOCCUS A
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified Ig - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.