

Product Datasheet

STREPTOCOCCUS PNEUMONIAE antibody (orb429823)

Description	Rabbit polyclonal antibody to STREPTOCOCCUS PNEUMONIAE
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Whole cell blend of numerous serotypes.
Target	STREPTOCOCCUS PNEUMONIAE
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	4.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.