

Product Datasheet

STREPTOCOCCUS PNEUMONIAE antibody (FITC) (orb429824)

Description	Rabbit polyclonal antibody against STREPTOCOCCUS PNEUMONIAE conjugated to FITC
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	Whole cell blend of numerous serotypes.
Target	STREPTOCOCCUS PNEUMONIAE
Preservatives	0.1% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	4 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.