

Product Datasheet

INFLUENZA A antibody (orb429827)

Description	Mouse monoclonal antibody to INFLUENZA A
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr
Immunogen	Influenza virus from Texas strain.
Target	INFLUENZA A
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.