

Product Datasheet

TROPONIN T CARDIAC ISOFORM 1 antibody (orb429828)

Description	Mouse monoclonal antibody to TROPONIN T CARDIAC ISOFORM 1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	HMW conjugate (FMPNLVPPKIC)
Target	TROPONIN T CARDIAC ISOFORM 1
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.