

Product Datasheet

BMP-9 antibody (orb429886)

Description	Rabbit polyclonal antibody to BMP-9
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	Synthetic peptide corresponding to the amino acids 87 - 98 (CKDDMGVPTLKYH) of the native molecule conjugated to keyhole limpet haemocyanin.
Target	BMP-9
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.