

Product Datasheet

CAMPYLOBACTER JEJUNI antibody (orb429903)

Description	Mouse monoclonal antibody to CAMPYLOBACTER JEJUNI
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr
Immunogen	Native, Type 1 Campylobacter jejuni.
Target	CAMPYLOBACTER JEJUNI
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgM - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.