

Product Datasheet

CANDIDA ALBICANS antibody (Biotin) (orb429914)

Description	Rabbit polyclonal antibody against CANDIDA ALBICANS conjugated to Biotin
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Biotin
Tested Applications	ELISA, ICC, IF, IHC-Fr
Immunogen	Candida albicans, type A
Target	CANDIDA ALBICANS
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	4.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	TYPE A
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.