

Product Datasheet

CHLAMYDIA TRACHOMATIS EBs antibody (orb429926)

Description	Goat polyclonal antibody to CHLAMYDIA TRACHOMATIS EBs
Species/Host	Goat
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Chlamydia trachomatis, strain L2, and other serovar groups.
Target	CHLAMYDIA TRACHOMATIS EBs
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.