

Product Datasheet

CHLAMYDIA TRACHOMATIS MOMP antibody (orb429927)

Description	Goat polyclonal antibody to CHLAMYDIA TRACHOMATIS MOMP
Species/Host	Goat
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr
Immunogen	Purified MOMP (major outer membrane protein) of Chlamydia trachomatis strain L2.
Target	CHLAMYDIA TRACHOMATIS MOMP
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified Ig - liquid
Concentration	4.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.