

Product Datasheet

CHLORAMPHENICOL antibody (orb429934)

Description	Sheep polyclonal antibody to CHLORAMPHENICOL
Species/Host	Sheep
Reactivity	Chemical
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Chloramphenicol base conjugated via its amino group to KLH using phenylenediisothiocyanate
Target	CHLORAMPHENICOL
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.