

Product Datasheet

hCG antibody (orb429944)

Description	Sheep polyclonal antibody to hCG
Species/Host	Sheep
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, RIA
Immunogen	Native, from urine.
Target	hCG
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	BETA
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.