

Product Datasheet

CHROMOGRANIN A antibody (orb429946)

Description	Sheep polyclonal antibody to CHROMOGRANIN A
Species/Host	Sheep
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	Synthetic peptide corresponding to the amino acids 1 - 15 (Ipnvspmknkgdtevm) of the native molecule conjugated to Keyhole limpet Haemocyanin.
Target	CHROMOGRANIN A
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.