

Product Datasheet

CHYMOTRYPSIN antibody (orb429947)

Description	Mouse monoclonal antibody to CHYMOTRYPSIN
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	Purified human pancreatic chymotrypsin.
Target	CHYMOTRYPSIN
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG3
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.