

Product Datasheet

CLOSTRIDIUM BOTULINUM TOXOID D antibody (orb429950)

Description	Goat polyclonal antibody to CLOSTRIDIUM BOTULINUM TOXOID D
Species/Host	Goat
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Clostridium botulinum D toxoid, formalin fixed.
Target	CLOSTRIDIUM BOTULINUM TOXOID D
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.