

Product Datasheet

COLLAGEN I antibody (orb429956)

Description	Rabbit polyclonal antibody to COLLAGEN I
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IH-Fr, IHC-P
Immunogen	Collagen type I from human placenta.
Target	COLLAGEN I
Preservatives	1% Mannitol1% Dextran. 0.1M Sodium Chloride0.01M Sodium Phosphate0.01M Sodium Borate
Form/Appearance	Purified IgG - lyophilized
Storage	-20°C after reconstitution
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.