

Product Datasheet

COLLAGEN I antibody (orb429965)

Description	Rabbit polyclonal antibody to COLLAGEN I
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P
Immunogen	Purified collagen type I from bovine dermis.
Target	COLLAGEN I
Preservatives	0.1% Sodium Azide (NaN ₃) Antibiotic antimycotic mixture 1%. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.