

Product Datasheet

COLLAGEN III antibody (orb429970)

Description	Rabbit polyclonal antibody to COLLAGEN III
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	Purified collagen III from rat skin.
Target	COLLAGEN III
Preservatives	0.1% Sodium Azide (NaN ₃) Antibiotic antimycotic mixture 1%. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.