

Product Datasheet

COLLAGEN I/II/III antibody (orb429971)

Description	Mouse monoclonal antibody to COLLAGEN I/II/III
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IH-Fr
Immunogen	Human placental collagen type 1.
Target	COLLAGEN I/II/III
Preservatives	1% Mannitol1% Dextran. 0.1M Sodium Chloride0.01M Sodium Phosphate0.01M Sodium Borate
Form/Appearance	Purified IgM - lyophilized
Storage	-20°C after reconstitution
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.